



STAY IN THE HOME  
**YOU LOVE**

[www.bruno.com](http://www.bruno.com)





## BRUNO'S BEST-SELLING STAIRLIFT, THE ELAN

Bruno's most popular stairlift can help restore full access to your home. Designed for straight indoor staircases, the Elan offers a mix of high-quality standard functionality, unsurpassed reliability, and ease-of-use design.

With Made-in-the-USA dependability, Bruno helps you stay in your home with comfort and ease.

Made-in-the-USA  
300 lb Lift Capacity (136 kg)

★★★★★  Verified order

"Bruno is an excellent, well-made product. We are pleased with our investment. The installation was easy and done professionally. American-made gives me confidence that it can take care of our needs for many years."

Ronald L.

## AFFORDABLE. QUIET. RELIABLE.

- **Simple Operation** – Control on armrest activates stairlift. Two wireless remote controls also included.
- **Space Saving** – Narrow vertical rail installs close to wall to preserve maximum stairway space.
- **Sleek & Clean** – Vertical rail with hidden gear rack gives clean appearance and reduces exposure to dirt, grease.
- **Easy On/Off** – Seat swivels up to 90 degrees at top landing, extends and locks in place for safe exit away from steps.
- **Compact** – Arms, seat and footrest fold to save space.
- **Dependable** – Battery-powered ensures performance even in power outage. Stairlift continually charging along entire rail.
- **Safe** – Easy to use retractable seat belt. Obstruction safety sensors stop lift when bumps into object.



Padded, generous-size seat.



Compact design for open space on stairs.



Verified order

“The Bruno stairlift has changed my life. It has made all parts of my home accessible to me. The chair is easy to operate, well constructed and safe. I love it!”

Kathleen C.



Offset swivel seat makes getting on/off easy.



Ergonomic armrest control.



Two wireless call/sends allow remote control access.



Sensors ensure safety.

## ELAN OPTIONS



Manual or power folding rail for narrow hallway or when doorway is at the bottom of stairs. Manual operation or push button automatic.



Power swivel seat for effortless exit. Armrest switch control or wireless call/send.



Power folding footrest automatically flips up/down when seat is raised/lowered. Arm switch control optional.



✓ Verified order

“Works smoothly and quietly to move my mother from floor to floor. Appreciate the safety features as well as the ease of operation.”

C. Savan



## TECHNICAL SPECIFICATIONS

**Seat Pad Size:** 15<sup>3/4</sup> in wide x 16<sup>3/4</sup> in deep  
(391 mm wide x 416 mm deep)

**Width Between Arms:** 20<sup>1/4</sup> in (514 mm)

**Footrest Size:** 10<sup>1/2</sup> in long x 11 in wide  
(267 mm long x 279 mm wide)

**Vertical Rail:**  
Installs within 5 in (127 mm) of the wall

## INDUSTRY-LEADING WARRANTY

**Lifetime Warranty on Major Components** -- motor, gearbox, drivetrain and rail -- for the lifetime of the Bruno Stairlift as long as the original purchaser owns the stairlift.

**2-Year Warranty** on all components after date of installation by an Authorized Bruno Dealer.

**30-Days Labor Coverage** after date of installation by an Authorized Bruno Dealer for any manufacturing defect.

## About Bruno

Family-owned Bruno Independent Living Aids has improved the lives of people with limited mobility for over 35 years. Engineered and manufactured in the USA, Bruno makes stairlifts and vertical platform lifts to help people better access homes, public buildings and businesses and scooter/powerchair lifts to assist people in easily transporting their mobility devices.

**Proudly made in the USA.**



Bruno Independent Living Aids, Inc.  
1780 Executive Drive, Oconomowoc, WI 53066  
[www.bruno.com](http://www.bruno.com)



CONSUMERAFFAIRS

© Bruno Independent Living Aids, Inc. 2019  
2019-08-03-Elan-6-PageV4.indd • MKT-SL-1