

## Aluminum Rehab Shower Commode Chair with Four Rear-locking Casters

HCPCS:



FIGURE A

- Versatile 3-in-1 design enables this model to easily be positioned over a standard toilet, used as a portable self-contained commode or shower chair.
- Simply and quickly adjust the height to achieve the ideal seated position.
- Padded flip-back armrests provide superior support and create a comfortable, non-slip surface for a safe and secure grip when sitting or standing.
- Easy to empty commode bucket to reduce splashes and spills. Use commode liners for even quicker clean up and less mess!
- 5" locking casters allow for extra security when the unit is not in motion.
- Durable and easy-to-clean padded seat (18 x 17) and backrest.
- Seat opening allows for unobstructed access to the perineal area for easy and effective cleansing.
- Swing-away detachable footrests with height and angle adjustments.
- Both lightweight and durable the rust-resistant aluminum frame supports up to 275 lbs.
- Simple, tool-free assembly.

### PRODUCTS

	Description
NRS185007	Shower&Commode Chair,1/ea

### SPECIFICATIONS

	Description
Back Height	19 "
Carton Shipping Weight	40 lbs.
Overall Height	39" - 41"
Overall Width Open	22.25"
Seat Dimensions	18" (W) x 17" (D)
Seat To Floor Height	20" - 22"
Warranty	Limited Lifetime
Weight	35.25 lbs.
Weight Capacity	275 lbs

### PARTS

	Description
STDS851-SA	Velcro Seat Belt Self Attaching,1ea, EA
NRS185006-21	BackLockingKnobForNRS185006/7, 1/pr, PR
NRS185006-14	DowntubePlasticInsrtsForNRS185006/7,1/ea, EA
NRS185006-12	FootplateForNRS185006/7, 1/pr, PR
NRS185006-16	FootrstPivotSaddleW/BoltNRS185006/7,1/ea, EA
NRS185006-22	AntiRattleInsertForNRS185006/7, 1/ea, EA

## PARTS

## Description

NRS185006-17	FootrestFrameLatchForNRS185006/7, 1/ea, EA
NRS185006-15	FootrstHgtAdjKnobW/BoltNRS185006/7,1/ea, EA
NRS185006-13	AngleAdjFtplateMntW/HdwrNRS185006/7,1/ea, EA
NRS185006-11	SwingawayFootrestForNRS185006/7, 1/pr, PR
NRS185006-09	FrameClipForNRS185006/7, 4/pkg, PKG
NRS185006-08	CasterWithLegForNRS185006/7, 1/ea, EA
NRS185006-05	ArmrestBracketLeftForNRS185006/7, 1/ea, EA
NRS185006-06	ArmrestBracketRightForNRS185006/7, 1/ea, EA
NRS185006-03	BackInsertForNRS185006/7, 1/ea, EA
NRS185006-04	ArmpadForNRS185006/7, 1/ea, PR
NRS185006-02	SeatClipW/ScrewsForNRS185006/7, 1/ea, EA
NRS185006-01	SeatWithInsertForNRS185006/7, 1/ea, EA