

VivaLIFT![®]
POWER RECLINER COLLECTION

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER

**NOW
AVAILABLE**
TREDELENBURG-LIKE
POSITIONING AND
MEMORY TOGGLE
HAND CONTROL!



Shown in Astro Brown Fabric

WITH
**COUNTLESS
POSITION
COMBINATIONS**

TRANQUIL 2

The Tranquil 2 combines quiet comfort with gorgeous fabrics, for a rich and upscale look. Get the support you need for your neck, shoulders and lower back as you dive into your favorite story.

The Tranquil 2 has been upgraded with an ergonomically designed memory toggle remote that features two programmable memory position buttons along with up/down/recline functions and custom comfort levels.

pridemobility.ca



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

TRANQUIL 2

VivaLIFT!
POWER RECLINER COLLECTION

COLLECTION	TRANQUIL 2	TRANQUIL 2	TRANQUIL 2	TRANQUIL 2
MODEL	PLR-935PW	PLR-935S	PLR-935M	PLR-935LT
POSITION	Infinite Relaxed Reclined Upward Lift Trendelenburg-like	Infinite Relaxed Reclined Upward Lift Trendelenburg-like	Infinite Relaxed Reclined Upward Lift Trendelenburg-like	Infinite Relaxed Reclined Upward Lift Trendelenburg-like
SUGGESTED HEIGHT*	5'3" and below	5'3" and below	5'4"-5'9"	5'10" and above
SIZE	Petite Wide	Small	Medium	Large/Tall
WEIGHT CAPACITY	400 lbs.	400 lbs.	400 lbs.	400 lbs.
SEAT-TO-FLOOR	18"	18"	19.5"	21.5"
SEAT DEPTH	18.5"	18.5"	20.5"	22.5"
SEAT WIDTH	24"	20"	22"	23.5"
TOP OF BACK TO SEAT	24.5"	25"	27"	30"
SITTING DEPTH**	36"	36.5"	38"	40"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort	Total Comfort	Total Comfort
CHAIR HEIGHT	38"	40"	42"	47"
DISTANCE FROM WALL	12.5"	13"	16"	16"
OVERALL WIDTH	36"	33.75"	37.5"	38"
CHAIR WEIGHT	143 lbs.	137 lbs.	145 lbs.	155 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated







* Suggested height range should be used only as a guide due to varying factors of each individual user.

** Seat depth from the front of the backrest to the end of the fully extended footrest.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



 = Indicates shortest available lead time

INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.
-  7. Trendelenburg-like position raises the seat with your feet elevated and your head reclined for the ultimate in relaxation and wellness.



STANDARD FEATURES

-  8. The memory button hand control adjusts all functions independently and allows you to charge your smart device.
-  9. Lithium battery backup ensures your recliner is operational when the power goes out.

STANDARD FABRIC: ASTRO



GREY

BROWN

MUSHROOM

OPTIONAL FABRIC: CRYPTON® ARIA



COOL GREY

ESPRESSO

Beamsville, ON / 888-570-1113 / pridemobility.ca

PRIDE
MOBILITY
LIVE YOUR BEST®