

# XPRESSO LITE CF Series



## EVOLUTION'S POPULAR XPRESSO LITE IN A CABLE-FREE VERSION.

With cable-free brakes, a center folding seat, durable folding basket, the Xpresso folds easily and is compact.

### FEATURES

**INNOVATIVE DESIGN** Roman arch provides support & stability

**SEAT** Soft, washable cushion

**BRAKE SYSTEM** Cable-free / Easy to adjust handle height

**COLLAPSIBLE BASKET** Foldable & removable

**WHEELS** 8" lightweight EVA tires

**BACK STRAP PAD** For added comfort and support

**COLORS** Ebony Black, Ruby Red, Azure Blue



### OPTIONAL ACCESSORIES

- Cup Holder
- Cane Holder
- Curb Climber
- Phone Holder
- Slow Down Brakes
- Eyeglass Holder
- Oxygen Tank Holder
- Wired Basket

### SPECIFICATIONS

#### MINI

15 lbs\*  
18" Seat Width  
18" Seat Height  
22.5" Width x 26" Length  
30.5" - 33.5" Handle Height  
350 lbs Weight Capacity

#### REGULAR

15 lbs\*  
18" Seat Width  
21" Seat Height  
22.5" Width x 26" Length  
30.5" - 34.5" Handle Height  
350 lbs Weight Capacity

#### TALL

15 lbs\*  
18" Seat Width  
24" Seat Height  
22.5" Width x 26" Length  
33.5" - 40" Handle Height  
350 lbs Weight Capacity

\*Net Weight of walker without accessories, backstrap and folding basket



Toll Free 1 800 556 2558

info@evolution.ca

www.evolution.ca