

## Airgo Fusion F20 Side-Folding Rollator & Transport Chair

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- The lightweight side-folding rollator and transport chair in one design
- Folds sideways easily by just a pull of the strap
- NEW covered soft basket folds with the frame and can be easily removed
- Patented Swing n' Lock footrests clip to the frame when folded
- The frame locks together when folded for safer lifting and storing
- Padded seat and a full backrest for superior comfort
- Dual armrests provide a place to rest yours arms in either mode
- Smooth bearings and an improved center of gravity for easy pushing
- Unique ergonomic handgrips fit naturally for easier brake access
- Features our new senior friendly tool-free brake adjustment knob
- Lightweight frame designed for superior stability
- BONUS curb climber included
- ADP# AA0000694
- Recognized by CAOT

PRODUCTS	Description
700-935	Airgo Fusion 2 in 1, 20", Blue, 1 EA
700-935BK	Airgo Fusion 2 in 1, 20", Black, 1 EA
700-935CR	Airgo Fusion 2 in 1, 20", Cranberry, 1EA

SPECIFICATIONS	Description
Handle Height	30" to 36.5"
Seat Height	20"
Seat Width	19.5"
Unit Width	27"
Warranty	Lifetime on Frame, 5 Years on Handle & Brakes, All other Parts 3 Months
Weight	17.1 lbs.
Weight Capacity	300 lbs

ACCESSORIES	Description
730-942	Slow brakes sfold, Black, 1PR, 25PR/CV, PR