

Chrome Knurled Grab Bar with Rotating Flange



HCPCS:




FIGURE A



FIGURE B

- Rotating flange for easy installation
- Knurled grab bar diameter provides a no-slip grip
- Mounts vertically or horizontally
- 3" distance from wall
- Supports up to 250lbs
- Limited Lifetime warranty

PRODUCTS

Description

12112RF	GrabBar,Chrm12"RotFlange,1/ea, RTL
12116RF	GrabBar,Chrm16"RotFlange,1/ea, RTL
12118RF	GrabBar,Chrm18"RotFlange,1/ea, RTL
12124RF	GrabBar,Chrm24"RotFlange,1/ea, RTL
12132RF	GrabBar,Chrm 32"RotFlange,1/ea, RTL