

Premium Plastic, Raised, Elongated Toilet Seat with Lock

HCPCS: E0244



FIGURE A



FIGURE B

- Designed for individuals who have difficulty sitting down or getting up from the toilet
- Adds 5" to height of toilet seat
- Tool-free arms are 21" apart and can be removed or added as needed
- Wide opening in front and back for personal hygiene
- Heavy-duty molded plastic construction provides additional strength and durability
- Standard locking mechanism ensures safety
- Lightweight and portable
- No tools required for installation
- Easy-to-clean
- Fits most elongated toilets
- Made using recycled materials

PRODUCTS

Description

12013	Raised Toilet Seat w/Arms,1/cs, RTL
-------	-------------------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	6.5 lbs.
Seat Dimensions	15.5" (W) x 17" (D) x 5" (H)
Warranty	Limited Lifetime
Weight	4.5 lbs.
Weight Capacity	300 lbs
Width Between Arms	21"