

Brochure

Salute RDX Mattress



Blake
Medical Distribution

1-866-662-3459
www.blakemedical.ca



The Salute RDX Mattress system provides alternating pressure therapy with micro low air loss to aid patients who are at risk of developing pressure ulcers or have existing stage 1, 2, 3, or 4 ulcers.

Caregiver friendly, the Salute RDX is an affordable, therapeutic surface solution for any clinical care setting. This micro low air loss mattress system utilizes 18 8" combined air cells (4.5" therapy and 3.5" static air) to ensure support even in the event of a power loss. The static air cell eliminates the need for a foam base, which means the whole surface can be easily disinfected.

Features

- 8" combined air cells
- Low air loss therapy
- 3-1 alteration
- Max inflate provides firm surface for transfers and wound dressing changes
- High vapor permeable and quilted nylon cover
- Happy Heel provides soft heel environment
- Integrated glide sheet to base cover for easy transferring
- Visual and audio alarm for low pressure and power failure with alarm mute function
- Keypad lock out function
- 8 Digital comfort setting
- Integrated power cable retainer
- SWL 500lbs
- Quick CPR deflation within 20 seconds



Happy Heel provides further therapy and comfort by decreasing pressure in patient's heel area



Control Unit

Size (LxWxH)	13.5" x 7.3" x 8.3"
Weight (lbs)	7.1
Cycle Time	5-20 min
Min Operating Pressure	12+/- 6mmHg
Max Operating Pressure	47 +/- 6mmHg
Flow-rate	8 liters per minute
Rated Voltage	AC 110-120V
Rated Frequency	60Hz
Protection Type	Class 1, Type BF
Max Current	0.2A
Low Pressure Alarm	Visual/Audio
Power Failure Alarm	Visual/Audio
Maintenance Service	LED Indicator
Functions	Static or 3-1 Alteration

Mattress

Size (LxWxH)	36" x 80" x 8"
Weight (lbs)	25.3
Cell Material	Blue Nylon with PU Coating
Cell Number	18 cells
Cover Material	Soft touch cover
Base Material	Black polyester with PVC Coating



