

# VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST  
COMFORTABLE  
RECLINER  
IN THE WORLD!

## THE 4-MOTOR RECLINER



FDA Class II Medical Device¹

Shown in Stonewash Granite Fabric



COUNTLESS  
POSITION  
COMBINATIONS

## URBANA 2¹

Innovative technology meets contemporary style for a recliner that supports you and your busy life. The **Urbana 2** by VivaLift!® Power Recliners is the perfect addition to your townhouse or condo.

The Urbana 2 has been upgraded with enhanced true-infinite positioning and an ergonomically designed, memory toggle remote that features two programmable memory position buttons, as well as up/down/recline functions and customized comfort.



Scan here for  
more information  
on VivaLift!®









**PRIDE**®  
MOBILITY  
LIVE YOUR BEST®

# URBANA 2<sup>1</sup>

**VivaLIFT!**<sup>®</sup>  
POWER RECLINER COLLECTION<sup>1</sup>

COLLECTION	URBANA 2 <sup>1</sup>
MODEL	PLR-965M
POSITION	True-Infinite Position
HEIGHT RANGE	5'10" and Above
SIZE	Medium
WEIGHT CAPACITY	400 lbs.
SEAT-TO-FLOOR	20"
SEAT DEPTH	20"
SEAT WIDTH	21"
TOP OF BACK TO SEAT	25"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	42"
DISTANCE FROM WALL	14"
OVERALL WIDTH	33"
CHAIR WEIGHT	145 lbs.
WARRANTY	7-Yr. Pro-rated


## INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. **NEW!** Enhanced true-infinite position recline enables the backrest and footrest to move independently so that a variety of positions can be reached including raising your feet above your heart.

## STANDARD FEATURES

- NEW!** Ergonomically designed toggle remote has two programmable memory position buttons. Performs the up/down, recline functions. Additionally, it has a convenient USB charger port that allows you to charge your smart device.
- Lithium battery backup ensures your recliner is operational when the power goes out.

STANDARD FABRIC: **STONEWASH** 



GRANITE

GUNMETAL

1. VivaLift!<sup>®</sup> Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



**ES** = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

(US) 800-800-8586 / (CAN) 888-570-1113 / [pridemobility.com](http://pridemobility.com)

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>