

Nitro Aluminum Rollator, 10" Casters

HCPCS:



FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Travel in style with the sleek design of the Euro-style Nitro Rollator. This rolling walker is cleverly engineered with a comfortable seat, convenient storage bag, adjustable height handles, and large 10 front casters for seamless rolling. With the brake cables inside the lightweight frame, the rollator walker provides a safe experience when you're out and about. With the ability to fold side-by-side, the Nitro Rollator is easily transported and stored.
- Adjustable contoured backrest
- Cross-brace design for added safety
- Clip hold and side-by-side folding
- Large and convenient storage bag included
- Weight capacity: 300 lb
- Weighs only 18 lb., easy-to-transport

PRODUCTS	Description
RTL10266	Nitro Rollator, Red Frame, 1/ea
RTL10266BK	Nitro Rollator, Black frame,1EA

SPECIFICATIONS	Description
Seat Dimensions	18" (W) x 10" (D)
Handle Height	33.5" - 38.25"
Weight Capacity	300 lb
Seat To Floor Height	20.5"
5 Year Limited Warranty	Brake and Handle Housing
Carton Shipping Weight	21.5 lb
Weight	18 lb

ACCESSORIES	Description
10266-CH	Nitro Cup Holder, 1/ea, EA
10266-BRAKE ENH	Nitro Slow Down Brake, ENH, 1PR, 25/CS, EA

PARTS	Description
1026602	HandGripLeft for Nitro,1/ea, EA
1026601-RN	HandBrakeAssembly,Left,10266,1ea, EA
1026603	HandGripRight for Nitro,1/ea, EA
1026601-W	HandBrakeAssembly,Left,10266WT,1/ea, EA
1026601-G	HandBrakeAssembly,Left,10266BK/10266BK-T, EA
1026630N	HandlebarForNitro/NitroTall,Side Mount, EA
1026604-B	BackrestAssemblyfor 10266 & 10266WT,1/ea, EA

PARTS	Description
1026605-W	BackrestLockingKnobW/Boltfor10266WT,1/pr, PR
1026605-G	BackrestLockingKnobW/Boltfor10266BK,1/pr, PR
1026604-G	Backrest Assembly for 10266BK,1/ea, EA
1026605-R	BackrestLockingKnobW/Boltfor10266,1/pr, PR
1026624	RightForkCompleteW/CasterforNitro, 1/ea, EA
1026625	LeftForkCompleteW/CasterforNitro, 1/ea, EA
1026617	FrontWheelWithBearingsforNitro, 1/ea, EA
1026616	WheelBrakeRightW/Hdwr,Nitro,Metal,1/ea, EA
1026615	WheelBrakeLeftW/Hdwr,Nitro,metal,1/ea, EA
1026612-R	RearWheelWithBearingsforNitro,Right,1/ea, EA
1026612-L	RearWheelWithBearingsforNitro,Left,1/ea, EA
1026609-FSW	SeatLiftingStrapW/AttachHrdwrefor10266WT, EA
1026609-FSG	SeatLiftingStrapW/AttachHrdwrefor10266BK, EA