

Muskoka LC-475

Home > Pride Power Lift Recliners > Muskoka LC-475




The **Muskoka LC-475** combines supportive comfort, smooth positioning, and practical everyday convenience. With enhanced three-position recline, power headrest, power lumbar, integrated USB A & C charging ports, and a fully padded chaise design, the LC-475 is built to help users relax, recharge, and rise with confidence.

Where to Buy

MSRP \$2,299 – \$2,499

Colour



Size

Small

MSRP \$2,299

Categories: Pride Power Lift Recliners, Sitting Pretty Power Lift Recliners

SKU: LC475S-DPH-H0A

Key Features

Where to Buy

Accessories

Additional Information

Comfort, Support, and Everyday Ease

The **Pride LC-475 Power Lift Chair** brings refined comfort and dependable lift assistance together in one thoughtfully designed recliner. Part of the **Muskoka Collection**, this chair is ideal for users who want a supportive

seating experience with premium comfort features built in.

Designed as an **enhanced three-position lift recliner**, the LC-475 makes it easy to sit, recline, relax, and transition to a standing position using a simple toggle hand control. Whether reading, watching TV, resting, or needing assistance getting up, this chair provides smooth, easy-to-use positioning for everyday comfort.

With a **fully padded chaise**, **adjustable waterfall back**, **power headrest**, and **power lumbar support**, the LC-475 is made to reduce pressure points while supporting the body from head to toe. Integrated charging ports, magnetic hand control docking, and battery backup add practical convenience for daily use.

LC-475

Muskoka Collection

COMFORT THAT LIFTS YOU

Designed for effortless support and personalised relaxation, the LC-475 brings comfort, independence, and style together in perfect balance.



MULTIPLE COMFORT POSITIONS

Find your perfect position



POWER HEADREST & LUMBAR

Custom support at the touch of a button



POWER LIFT ASSIST

Smooth, steady lift for safe standing



PREMIUM FABRIC

Soft, durable, and easy to maintain



EASY-TO-USE REMOTE
Simple controls for headrest, lumbar, recline, and lift functions.

MULTIPLE COMFORT POSITIONS

One chair. Every position you need.



RELAX
Kick up your feet and unwind.



RECLINE
Lie back and enjoy full-body comfort.



LIFT
Gentle lift assistance for confidence and independence.



STAND
Smooth transition to help you stand with ease.

Key Features

Enhanced Three-Position Recline

The LC-475 offers enhanced three-position functionality, giving users comfortable positioning options for sitting upright, relaxing, or reclining. The lift function helps support a smoother transition from seated to standing.

Easy-to-Use Toggle Hand Control

The simple toggle hand control allows users to adjust the chair with ease. Its intuitive design makes it straightforward to operate the lift and recline functions without complicated buttons or settings.

Power Headrest

The built-in power headrest allows users to adjust head and neck positioning for improved comfort while reading, watching TV, or relaxing. This is especially useful for maintaining better line of sight while reclined.

Power Lumbar Support

Power lumbar adjustment provides added lower-back support, helping users fine-tune the chair to their comfort preferences throughout the day.



Specifications

Specification	LC-475S Small	LC-475M Medium	LC-475L Large
Collection	Muskoka Collection	Muskoka Collection	Muskoka Collection
Position	Enhanced 3-Position	Enhanced 3-Position	Enhanced 3-Position
Suggested Height Range	5'3" & below	5'4" – 5'9"	5'9" – 6'
Weight Capacity	400 lb	400 lb	400 lb
Back Style	Waterfall, adjustable	Waterfall, adjustable	Waterfall, adjustable
Seating Type	Total Comfort	Total Comfort	Total Comfort
Footrest Extension	Standard	Standard	Standard
Fully Padded Chaise	Yes	Yes	Yes
Seat-to-Floor Height	18.5"	20"	22"
Seat Width	20"	22"	24"
Seat Depth	18"	20"	22"
Seat Depth with Footrest Extended	39"	41"	43.5"
Top of Back to Seat	24"	26"	27"
Overall Width	33"	34.5"	35.5"

Fully Padded Chaise Design

The chaise pad design fully supports the legs and helps reduce pressure points. This creates a more continuous, supportive seating surface from the seat through the footrest.

Adjustable Waterfall Back

The LC-475 features a comfortable waterfall-style back that is adjustable, allowing users to enjoy plush, contoured support.

Integrated USB A & C Charging Ports

Convenient USB A and USB C charging ports make it easy to keep phones, tablets, and other compatible devices powered while relaxing.

Magnetic Hand Control Docking on Both Sides

The hand control can be docked magnetically on either side of the chair, making it easier to keep the controller within reach and reducing the chance of it slipping between cushions or falling to the floor.

Integrated Battery Backup

The integrated battery backup provides added peace of mind in case of a power failure, helping ensure the chair can still be operated when needed.

Available in Multiple Sizes

The LC-475 is available in **Small, Medium, and Large** sizes, making it easier to choose a chair that fits the user's height, comfort needs, and room layout.